

# REDUCE YOUR WATER FOOTPRINT

PATENTED ENVIROVAC® WATER CONSERVATION

The patented EnviroVac® system is a patented water conservation system engineered for optional integration into all Beta Star sterilizers. This system reduces the amount of water used during a sterilization cycle by as much as 75%. Based on cycle per day usage rates, you can easily recover your cost within the first six months.

The performance of the EnviroVac® is not affected by normal fluctuations in water temperature or supply pressure. This ensures you the most reliable and environmentally friendly engineered sterilizer available.

## ENVIROVAC® BENEFITS

- Reduces sterilizer water consumption up to 75%
- Lowers operating costs
- Reduces wastewater
- Conserves energy and water
- Used for green building upgrades
- Demonstrates environmental stewardship
- Reduces water treatment costs
- Gains tax benefits for the facility
- Achieves LEED ratings

20" x 20" x 38" Load Size	EnviroVac® Total Gallons	Traditional Pump Total Gallons
Small	56	226
Medium	64	256
Large	70	286

24" x 36" x 48" Load Size	EnviroVac® Total Gallons	Traditional Pump Total Gallons
Small	126	330
Medium	142	346
Large	154	358

**BETA STAR**  
LIFE SCIENCE EQUIPMENT

