

# OIL FREE AIR COMPRESSOR



## THE EQUIPMENT


The oil free air compressor, available from Beta Star, provides an alternative to house compressed air for existing facilities. The oil free design requires minimal service to keep the compressor running smoothly.

## FEATURES

- Twin domed head assemblies
- Compressor body made entirely out of metal
- Cabinet installation ready
- 50µm air intake filter
- 5µm filter regulator with automatic drain
- Corrosion resistant powder coated receiver
- Biased valves and unloader valve for easy restart
- Twin seals allow higher pressure with less noise

## COMPRESSOR SPECIFICATIONS

Specification	Label	Number
Voltage	V	1x120
Frequency	Hz	60
Motor	HP	1/2
Displacement	CFM	3.2
Air flow @ 120 PSI	CFM	1.2
Max. pressure	PSI	120
Max. current	A	5.7
Receiver volume	Gallons	1.1
Noise level @ 1m	dB(A)	63
Weight	lbs	40
Dimensions	in	15.9 x 12.0 x 13.3



**DECIBEL RANGE: 48 - 75**

