

STERILIZER LOADING EQUIPMENT



LOADING CART AND TRANSFER CARRIAGE

The 316L stainless steel loading cart is designed to hold the goods or media being sterilized. For smaller sterilizers, the loading cart rides on the tracks between the transfer carriage and the vessel. The transfer carriage is used to transport the loading cart from station to station. A drip pan can be added to the loading cart to capture spilled liquid or media.

For pit mounted bulk (vivarium) sterilizers, large loading carts holding hundreds of cages can be rolled directly into the sterilization chamber. While cage throughput is defined by your sterilization chamber and cage size, our sales team can provide an accurate cage per cycle number. To get your customized cage per cycle throughput, call us today at **610.273.4222** or visit www.betastar.com.



24" X 36" LOADING CART AND TRANSFER CARRIAGE

SMALL TO BULK LOADING OPTIONS

